



CLASSIFICATION OF WHEELCHAIR FOR CLAMPING RINGS

BRAND	MODEL	TUBE SHAPE	FRAME DIAMETER	TUBE FRAME	CLAMPING RINGS
4ma3ma	Tuk	Round	25,4	parallel / convergent	AGP25 / AGCP125
Aria	1.0	Round	30,4	parallel	AGP32 + RID32-30-0000
	Speciale	Round	30,4	parallel	AGP32 + RID32-30-0000
	Ultra	Round	30,4	parallel	AGP32 + RID32-30-0000
	2.0	Round	30,4	parallel	AGP32 + RID32-30-0000
Drive Medical	Rotec	Round	25,4	parallel	AGP25
Etac	Cross 5 / 5 XL	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Cross 6	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Crissy Active	Round	**to check	parallel	**to check
	Crissy Swing-away (*)	Round	25,4	parallel	AGP25
	(*) the Clamping Rings have to be installed on the bottom tubes (where are fixed the front wheel - see Installation Wheelchair Catalogue).				
Excel	G-Lite Pro	Round	22,5	parallel	AGP25 + RID25-22-0000
	(*) the Clamping Rings have to be installed on the bottom tubes (where are fixed the front wheel - see Installation Wheelchair Catalogue).				
Jumper	WCMX Parkour	Round	32,4	parallel	AGP32
	WCMX Pro Action	Round	32,4	parallel	AGP32
Ki Mobility	Catalyst 5	Round	25,4	parallel	AGP25
	Catalyst 5Ti	Round	25,4	parallel	AGP25
	Catalyst 5Vx	Round	25,4	parallel	AGP25
	Catalyst 4 / 4C (*)	Round	28,6	parallel	AGP32 + RID32-28-0000
	Rogue	Round	34,9	parallel	AGP35
	Rogue ALX	Round	34,9	parallel	AGP35
	Tsunami ALX (*to check availability)	Round	34,9	parallel	AGP35
	Rogue XP	Round	34,9	parallel	AGP35
	Ethos	Round	28,6	parallel	AGP32 + RID32-28-0000
	Little Wave Arc	Round	28,6	parallel	NOT AVAILABLE
	Little Wave Flip	Round	28,6	parallel	NOT AVAILABLE
	Little Wave Cliek	Round	28,6	parallel	AGP32 + RID32-28-0000
	Little Wave XP	Round	28,6	parallel	AGP32 + RID32-28-0000

	Spark	Round	28,6	parallel	AGP32 + RID32-28-0000
(*) the Clamping Rings have to be installed on the bottom tubes (where are fixed the front wheel - see Installation Wheelchair Catalogue).					

Küschall	Compact Attract 1.0	Round	30,4	parallel	AGP32 + RID32-30-0000
	K-Series Attract 1.0	Round	30,4	parallel	AGP32 + RID32-30-0000
	KSL 1.0 / 2.0 (*)	Round	30,4	parallel	AGP32 + RID32-30-0000
	Champion 1.0 / 2.0 (*)	Round	30,4	parallel	AGP32 + RID32-30-0000
	Compact 2.0	Round	25,4	parallel	AGP25
	K-Series 2.0	Round	30,4	parallel	AGP32 + RID32-30-0000
	Advance (until end of 2019)	Round	30,4	parallel	AGP32 + RID32-30-0000
	R-Series (until end of 2019)	Round	30,4	divergent	AGCP132 + RID32-30-0000
	Compact 1.0 (until May 2020)	Round	30,4	parallel	AGP32 + RID32-30-0000
	Ultra Light 1.0 (until May 2020)	Round	30,4	parallel	AGP32 + RID32-30-0000
	K-Series 1.0 (until May 2020)	Round	30,4	parallel	AGP32 / AGCP132 + RID32-30-0000
	K-Series Titanium (until May 2020)	Round	28,4	parallel	AGP32 + RID32-28-0000
	K-Series Carbon (until May 2020)	Round	30,4	parallel	NOT AVAILABLE
Champion Carbon 1.0 (until May 2020)	Round	30,4	parallel	NOT AVAILABLE	
(*) In order to install the KLICK's Clamping Rings, the wheelchair in question should be adapted by cutting the seat around the clamping rings (see Installation Wheelchair Catalogue).					

Icon	A1 (*)	Round	25,4	parallel	AGP25
(*) the Clamping Rings have to be installed on the bottom tubes (see Installation Wheelchair Catalogue).					

Mediland	Kometa	Round	25,4	parallel	AGP25
----------	--------	-------	------	----------	-------

Meyra	Hurricane Daily	Round	25,4	parallel / convergent	AGP25 / AGCP125
	Hurricane S	Round	25,4	parallel / convergent	AGP25 / AGCP125
	Fuse R	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Nano S	Round	25,4	parallel	AGP25
	Nano C	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Nano X	Round	25,4	parallel	AGP25
	Nano	Round	25,4	parallel / convergent	AGP25 / AGCP125
	Smart S	Round	25,4	parallel	AGP25
	Avanti (*)	Round	25,4	parallel	AGP25
	Avanti Pro	Round	25,4	parallel	AGP25
	Smart F (*)	Round	25,4	parallel / convergent	AGP25 / AGCP125
(*) the Clamping Rings have to be installed on the bottom tubes (see Installation Wheelchair Catalogue).					

Molab	G6	Round	30,4	parallel	AGP32 + RID32-30-0000
	G6 Kid	Round	25,4	parallel / divergent	AGP25 / AGCP125
Motion Composite	Apex A	Round	32,4	parallel	AGP32
	Apex C	Round	32,4	parallel	AGP32
	Helio C2	Round	32,4	parallel	AGP32
	Helio C2 Carbon	Round	32,4	parallel	AGP32
	Veloce Carbon	Round	32,4	parallel	AGP32
	Move	Round	32,4	parallel	AGP32
	Helio A6	Round	32,4	parallel	AGP32
	Helio A7	Round	32,4	parallel	AGP32
Nuova Blandino	Torino	Round	25,4	parallel	AGP25
	GR120	Round	25,4	parallel	AGP25
Offcarr	Althea	Round	25,4	parallel	AGP25
	Alhena	Round	25,4	parallel	AGP25
	Halley	Round	25,4	parallel	AGP25
	Diva	Round	25,4	parallel	AGP25
	Eos	Round	31,4	parallel	AGP32
	Eos 3	Round	31,4	convergent	AGCP132
	Idra 2.0	Round	25,4	parallel	AGP25
	XLR8	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Fixed	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Fenice	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Themis	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Venus	Round	31,4	parallel	AGP32
	Venus Adventure	Round	31,4	parallel	AGP32
	Venus Elite	Round	31,4	parallel	AGP32
	Quasar / Quasar GT	Round	25,4	parallel	AGP25
	QPX	Round	25,4	parallel	AGP25
	Zodiac	Round	31,4	parallel	AGP32
	Vega	Round	25,4	parallel	AGP25
Vega E	Round	25,4	parallel	AGP25	
	Avantgarde 4 DV, DS, CLT	Round	25,4	parallel	AGP25
	Avantgarde CV, CS, CLT	Round	25,4	parallel	AGP25

Ottobock	Start4 M1 / M2 (all models)	Round	25,4	parallel	AGP25
	Motus CS / CV / 2.0	Round	25,4	parallel	AGP25
	Avantgarde XXL 2	Round	25,4	parallel	AGP25
	Invader (*)	Round	25 / 35	parallel / convergent	AGP25 / AGCP135
	Invader Wave	Round	30	parallel	NOT AVAILABLE
	Zenit & Zenit R (Carbon)	Round	30	parallel	AGP32 + RID32-30-0000
	Zenit & Zenit R (Aluminium)	Round	30	parallel	AGP32 + RID32-30-0000
	Zenit CLT & Zenit R CLT	Round	30	parallel	AGP32 + RID32-30-0000
	Ventus	Round	30	parallel	AGP32 + RID32-30-0000
	Voyager	Round	35	parallel	AGP35
	Avantgarde Teen2/VR	Round	Custom	parallel	NOT AVAILABLE
	Bravo Racer	Round	Custom	parallel	NOT AVAILABLE
	Blizzard	Round	35	parallel	NOT AVAILABLE
(*) the Clamping Rings AGP25 are for the installation on the seat tube, the clamping rings AGCP135 are for the installation on the lower tubes.					

Panthera	S3 Swing Familjen (*)	Round	22	convergent	AGCP125 + RID25-22-0000
	S3 Familjen (*)	Round	22	convergent	AGCP125 + RID25-22-0000
	U3 (*)	Round	22	convergent	AGCP125 + RID25-22-0000
	U3 Light (*)	Round	22	convergent	AGCP125 + RID25-22-0000
	X (Carbon)	Oval	OVAL Custom	parallel	NOT AVAILABLE
(*) the Clamping Rings have to be installed on the bottom tubes (see Installation Wheelchair Catalogue).					

Pro Activ	Speedy 4teen (*)	Round	25,4	convergent	AGCP125
	Litty (all models) (*)	Round	25,4	parallel	AGP25
	Buddy (*)	Round	25,4	parallel	AGP25
	Speedy (all models)	to check			
	Traveler (all models)	Round	25,4	parallel	AGP25
(*) the Clamping Rings have to be installed on the bottom tubes (see Installation Wheelchair Catalogue).					

Progeo	Exelle vario	Round	25,4	parallel	AGP25
	Tekna Advance	Oval	OVAL 25x35	parallel	AGPOV
	Tekna Advance Carbon	Oval	OVAL 25x35	parallel	AGPOV
	Yoga	Oval	OVAL 25x35	parallel	AGPOV
	Ego/ Ego Custom	Oval	OVAL 25x35	parallel	NOT AVAILABLE
	Ego Carbon	Oval	OVAL 25x35	parallel	NOT AVAILABLE
	Joker / Joker Energy	Oval	OVAL 25x35	parallel	AGPOV
	Joker Full Carbon	Oval	OVAL 25x35	parallel	AGPOV

	Joker Carbon Pack	Oval	OVAL 25x35	parallel	AGPOV
	Joker Evolution	Oval	OVAL 25x35	parallel	AGPOV
	Joker R2	Oval	OVAL 25x35	parallel	AGPOV
	Carbomax	Oval	OVAL Custom	parallel	NOT AVAILABLE
	Duke	Oval	OVAL Custom	parallel	AGPOVPC
	Noir 1.0	Oval	OVAL Custom	parallel	AGPOV
	Noir 2.0	Oval	OVAL Custom	parallel	AGPOVPC

RGK	Tiga (all models)	Round	25,4	parallel / convergent	AGP25 / AGCP125
	Veypr Sub4	Oval	OVAL Custom	parallel	NOT AVAILABLE
	Octane (all models)	Round	32,4	parallel / convergent	AGP32 / AGCP132
	Maxlite	Round	19,0	convergent	NOT AVAILABLE
	Hi Lite (all models)	Round	25,4	parallel / convergent	AGP25 / AGCP125

Schmicking	F-Design (Foldable)	Round	25,4	parallel / convergent	AGP25 / AGCP125
	L-Design	Round	35	parallel / convergent	AGP35
	X-Design	Round	25,4	parallel / convergent	AGP25 / AGCP125

Sorg	Mio	Round	30,4	parallel / convergent	AGP32 / AGCP132 + RID32-30-0000
	Vector /Vector BSA	Round	25,4	parallel	NOT AVAILABLE
	Jump (all models)	Round	25,4	parallel	NOT AVAILABLE
	Airon	Round	30,4	parallel	AGP32 + RID32-30-0000

Sunrise (Quickie)	Helium	Oval	OVAL 27x33	parallel	AGPOVSX
	Helium Freestyle	Oval	OVAL 27x33	parallel	AGPOVSX
	Nitrum	Oval	OVAL 27x33	parallel	AGPOVSX
	Krypton R	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Krypton F	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	Argon (all models)	Round	28	parallel	AGP32 + RID32-28-0000
	Neon	Round	28	parallel	AGP32 + RID32-28-0000
	Neon 2	Round	28	parallel	AGP32 + RID32-28-0000
	Life	Round	23	parallel	AGP25 + RID25-23-0000
	Life R	Round	23	parallel	AGP25 + RID25-23-0000
	Life T	Round	23	parallel	AGP25 + RID25-23-0000
	Helix 2	Round	25,4	parallel	AGP25
	QS5 X	Round	25,4	parallel	AGP25
	Xenon (all models)	Oval	OVAL 27x33	parallel	AGPOVSX

	Easy (all models)	Round	23	parallel	AGP25 + RID25-23-0000
TiLite	ZR (all series)	Round	31,8	parallel / convergent	AGP32 / AGCP132
	ZRA (all series)	Round	31,8	parallel / convergent	AGP32 / AGCP132
	TR (all series)	Round	25,4	parallel / convergent	AGP25 / AGCP125
	TRA	Round	25,4	parallel / convergent	AGP25 / AGCP125
	2GX (all series)	Round	25,4	parallel	AGP25
	TX (until 2017)	Round	25,4	parallel	AGP25
	AEROZ (all series)	Round	31,8	parallel / convergent	AGP32 / AGCP132
	AEROT	Round	31,8	parallel / convergent	AGP32 / AGCP132
	AEROX (all series)	Round	25,4	parallel	AGP25
	TWIST	Round	25,4	parallel / divergent	AGP25 / AGCP125
Vassilli	17.68N	Round	25,4	parallel	AGP25
	17.70	Round	25,4	parallel	AGP25
Vermeiren	Trigo	Oval	Custom	parallel	AGPOVT
	D200-V	Round	25,4	parallel	AGP25
	V300D	Round	25,4	parallel	AGP25
	Bobby / Bobby Evo	Round	25,4	parallel	NOT AVAILABLE
	D100	Round	25,4	parallel	AGP25
Wolturnus	Hawk (*)	Round	25,4	parallel	AGP25
	W5 (*)	Round	25,4	parallel	AGP25
	W5 SL (*)	Round	22,4	parallel	AGP25 + RID25-22-0000
	W5 S (*)	Round	25,4	parallel	AGP25
	W5 XXL (*)	Round	25,4	parallel	AGP25
	W5 K (*)	Round	25,4	parallel	AGP25
	W5 D (*)	Round	25,4	parallel	AGP25
	Tukan (*)	Round	30,4	parallel	AGP32 + RID32-30-0000
	Tukan A (*)	Round	30,4	parallel	AGP32 + RID32-30-0000
	A-Run	Round	25,4	parallel	NOT AVAILABLE
	Merlin (foldable)	Round	25,4	parallel	AGP25
	Gitano	Custom	NOT AVAILABLE	parallel	NOT AVAILABLE
	(*) In order to install the KLiCK's Clamping Rings, the straps fo the seats have to be positionated around the clamps (see Installation Wheelchair Catalogue).				